


ACCEL[®] All-Purpose and
Biliary Drainage Catheter Sets
with TrueGlide[®] Hydrophilic Coating


ACCEL[®] Secure
Catheter Stabilization Device

Designed and Developed for Optimal Patient Comfort and Clinical Care

The ACCEL® All-Purpose and Biliary Drainage Catheters are designed to maximize fluid drainage volume with five large, oval drainage holes and deliver a smooth atraumatic catheter insertion with the TrueGlide® Hydrophilic Coating.


Our comprehensive portfolio of sizes and configurations is designed to meet all your patients' drainage needs.





Clinical Features:

- Used for biliary, abscess and other multipurpose drainage
- Available in sizes 5-8 Fr, 10 Fr, 12 Fr, 14 Fr and 16 Fr with French size printed on strain relief
- Available in catheter working lengths 15 cm, 20 cm, 25 cm, 35 cm and 40 cm
- Allows for 0.035 in. (5-7 Fr) and 0.038 in. (8-16 Fr) guidewire diameter
- All-purpose drainage catheter placement can be performed using direct access or the Seldinger Technique
- Flexible polyurethane material and reinforced tip designed to help prevent kinking and provide optimal catheter placement



ACCEL® Closed Loop Locking Pigtail All-Purpose Drainage Catheter Set

Set contents:

- One - ACCEL All-Purpose Drainage Catheter
- One - Splittable Straightener
- One - Trocar Stylet
- One - Metal Stiffening Cannula
- One - Flexible Stiffening Cannula (not included with 6 Fr set)
- One - Dead End Cap

Reference Number	French Size	Working Length (cm)	Straightened Catheter Length (cm)	Radiopaque Marker Band	Pigtail Drain Holes	Case Qty.
613901	6 Fr*	25	30		5	5
613902	8 Fr	25	30		5	5
613903	10 Fr	25	30		5	5
613904	12 Fr	25	30		5	5
613905	14 Fr	25	30		5	5
613906	16 Fr	25	30		5	5
613907	8 Fr	25	30	✓	5	5
613908	10 Fr	25	30	✓	5	5
613909	12 Fr	25	30	✓	5	5
613910	14 Fr	25	30	✓	5	5
613911	8 Fr	35	40		5	5
613912	10 Fr	35	40		5	5
613913	12 Fr	35	40		5	5
613914	14 Fr	35	40		5	5
613915	8 Fr	40	45	✓	5	5
613916	10 Fr	40	45	✓	5	5
613917	12 Fr	40	45	✓	5	5
613918	14 Fr	40	45	✓	5	5

*0.035 in. guidewire compatibility
0.038 in. guidewire compatibility

Available with or without
Radiopaque Marker Band



Excellent Imaging Properties
with and without radiopaque
marker band under fluoroscopy

ACCEL® MINI-Closed Loop Pigtail All-Purpose Drainage Catheter Set

Set contents:


- One - ACCEL MINI All-Purpose Drainage Catheter
- One - Splittable Straightener
- One - Trocar Stylet
- One - Metal Stiffening Cannula
- One - Flexible Stiffening Cannula (not included with 5 Fr and 6 Fr sets)
- One - Dead End Cap



Reference Number	French Size	Working Length (cm)	Straightened Catheter Length (cm)	Locking Mechanism	Pigtail Drain Holes	Case Qty.
613919	5 Fr*	15	20		4	5
613920	6 Fr*	15	20		4	5
613921	7 Fr*	15	20		4	5
613922	6 Fr*	15	20	✓	4	5
613923	8 Fr	15	20	✓	4	5
613924	10 Fr	15	20	✓	4	5
613925	12 Fr	15	20	✓	4	5
613926	7 Fr*	25	30	✓	4	5
613927	8 Fr	25	30	✓	4	5
613928	10 Fr	25	30	✓	4	5
613929	12 Fr	25	30	✓	4	5
613930	14 Fr	25	30	✓	4	5

*0.035 in. guidewire compatibility
0.038 in. guidewire compatibility

Locking and Non-locking
Configurations



ACCEL® Closed Loop Locking Pigtail Biliary Drainage Catheter Set with Radiopaque Marker Band

Set contents:

- One - ACCEL Biliary Drainage Catheter
- One - Splittable Straightener
- One - Metal Stiffening Cannula
- One - Flexible Stiffening Cannula
- One - Dead End Cap

Reference Number	French Size	Working Length (cm)	Straightened Catheter Length (cm)	Radiopaque Marker Band	Pigtail Drain Holes	Case Qty.
613931	8 Fr	40	45	✓	5	5
613932	10 Fr	40	45	✓	5	5
613933	12 Fr	40	45	✓	5	5
613934	14 Fr	40	45	✓	5	5

0.038 in. guidewire compatibility



ACCEL® Universal J-Curved Non-Locking Tip All-Purpose Drainage Catheter Set

Set contents:


- One - ACCEL J-Curved All-Purpose Drainage Catheter
- One - Splittable Straightener
- One - Trocar Stylet
- One - Metal Stiffening Cannula
- One - Flexible Stiffening Cannula (not included with 6 Fr set)
- One - Dead End Cap

Reference Number	French Size	Working Length (cm)	Straightened Catheter Length (cm)	J-Curve Drain Holes	Case Qty.
613935	6 Fr*	20	25	5	5
613936	8 Fr	25	30	5	5
613937	10 Fr	25	30	5	5
613938	12 Fr	25	30	5	5

*0.035 in. guidewire compatibility
0.038 in. guidewire compatibility

ACCEL Secure Catheter Stabilization Device

ACCEL Secure is breathable and enhances patient comfort while maintaining catheter fixation and minimizing site leakage. Hydrocolloid adhesive with an outer adhesive cover absorbs superfluous exudates and eases device removal.



Reference Number	French Size	Case Qty.
613942	5-12 Fr	25
613943	12-18 Fr	25

Indications for Use:

ACCEL® All-Purpose Drainage Catheter Sets with TrueGlide® Hydrophilic Coating are designed for percutaneous drainage of abscesses applied in the gastroenterological system, urological system and other applicable drainage abscesses.

ACCEL® Biliary Drainage Catheter Sets with TrueGlide Hydrophilic Coating are designed for percutaneous drainage of abscess fluid, cyst, gallbladders, nephrostomy, urinary and other fluids.

ACCEL® Secure Catheter Stabilization Device is a sterile device for percutaneous catheters and tubes, consisting of a hydrocolloid adhesive with an outer adhesive cover. Refer to Instructions for Use for relevant warnings, precautions, complications and contraindications.

For more information or to place an order, contact your B. Braun Interventional Systems Inc. representative or call 877-836-2228

Rx only

ACCEL and TrueGlide are registered trademarks of B. Braun Interventional Systems Inc.

©2024 B. Braun Interventional Systems Inc. BIS16_3/4/24



Not made with natural rubber latex



Not made with DEHP



Not made with PVC

Distributed by:

B. Braun Interventional Systems Inc. | Part of the B. Braun Group of Companies

Bethlehem, PA 18018 | USA | Tel 877-836-2228 | Fax 610-849-1334 | www.bisusa.com